













ALVO Rapsodia

Perfect for following disciplines:

-  Bariatric Surgery (BS)
-  Cardiovascular (CV)
-  Colorectal Surgery (CRS)
-  General Surgery (GS)
-  Gynecology (GYN)
-  Minimally Invasive Surgery (MIS)
-  Neurosurgery (NS)
-  Ophthalmology (OPHT)
-  Orthopedic and Traumatology (OT)
-  Otorhinolaryngology (ENT)
-  Spinal Surgery (SS)
-  Urology (UROL)



* Available options:

- powered leg support
- alternative drive

Remote controlled functions:



Height



Trendelenburg /
Reverse
Trendelenburg



Lateral
tilt



Longitudinal
shift



Back
section



Leg
support *



Tabletop
leveling



Flex
position



Reflex
position



Brake

Universal mobile surgical table with highly competitive specifications

Working ergonomics

- ergonomic base design to provide unrestricted surgical access with exchangeable leg section and slim base at the side
- smooth edges to stop liquids from dropping on the floor




Ease of configuration

- quick set-up and segments exchange thanks to click-in connection of segments




Intraoperative pressure ulcers prevention

- with waterproof ultrasonic-welded soft pads



Surgical Site Infections prevention

- with Velcro-free table pads and sleek side rails for easier cleaning



Intuitive control interfaces

- intuitive control with backlit remote controller (both standard and wireless) or alternatively via OR integration control panel (ALVO Integra)



Universal mobile surgical table with highly competitive specifications



Wide adjustment
ranges for various
surgical disciplines

